


OFFERTA PROMOZIONALE

Cytoskeleton Signal Seekers 15% off

La Società Italiana Chimici, in qualità di Distributore della Cytoskeleton Inc. offre uno sconto del 15% su tutti i prodotti presenti sul flyer:



L'offerta promozionale è relativa alle seguenti categorie di prodotti:

- **Signal Seeker™ Ubiquitin Tools**
- **Signal Seeker™ SUMOylation 2/3 Tools**
- **Signal Seeker™ Phosphotyrosine Tools**

CONDIZIONI DI VENDITA:

Prezzi: al netto dell'IVA (22%)

Spese di Trasporto: € 25,00 (per ordini inferiori ad € 500,00)

Consegna: 10 giorni lavorativi dal ricevimento dell'ordine

Validità dell'offerta: fino al 31 Ottobre 2026

Codice Offerta: cyto#4-2026

CONTATTARE:

Società Italiana Chimici

tel: +39 06 8818936 - 8800211

e-mail: sales@sichim.com



Signal-Seeker PTM Detection Kits



Signal-Seeker Toolkits

- Comprehensive kits
- Minimal optimization time
- Maximizes PTM detection
- New Acetyl Lysine Affinity Beads and Kits

 **Keep It Simple**
Signal-Seeker™ PTM Toolkits

[Click to learn more](#)

15% Discount on Select Signal-Seeker Products



Product Summary Highlights


Acetyl-Lysine antibody (Cat. # AAC04-Beads)

Simplified acetylation detection with Cytoskeleton's Acetyl Lysine Affinity Beads-AAC04

Fig Legend: AAC03 & AAC02 acetyl lysine affinity beads (50µl bead slurry) were used to IP acetylated proteins from Cos-7 cells either treated (+) or untreated (-) with deacetylase inhibitors.

Learn more at: <https://www.cytoskeleton.com/product/acetyl-lysine-affinity-beads-aac04>

Measuring total acetylation levels in cell and tissue extracts




SUMO1 affinity beads (Cat. # ASM11-beads)

The only SUMO1 affinity beads commercially available with precise specificity for SUMO1 modified proteins

Investigate endogenous SUMO1 modifications for any target protein

Learn more at: <https://www.cytoskeleton.com/product/sumoylation-1-affinity-beads>



Signal-Seeker Ubiquitin Detection Kit (Cat. # BK161)

The original, comprehensive ubiquitin detection kit: contains all of the essential components needed for effective ubiquitination detection

Perfect for scientists that are new to the Ubiquitin field

Learn more at: <https://www.cytoskeleton.com/product/signal-seeker-ubiquitin-detection-kit>

Essential Kit Components:

- IP Affinity Beads: (Cat. # UBA01-beads)
- IP control Beads: (Cat. # CUB02-beads)
- Universal Lysis System
- Key Inhibitors
- Wash Reagents
- Elution Reagents
- Spin Columns
- Protein Quantitation System
- Chemiluminescence

Cytoskeleton Inc. Campaign #4 2026 Signal Seekers

CAT. NO.	PRODUCT	AMOUNT	Shipping Temp.	Storage on Arrival
AAC01-L	Acetyl Lysine Antibody Mouse Monoclonal	2 x 100 µl	AT	4°C
AAC01-S	Acetyl Lysine Antibody Mouse Monoclonal	1 x 25 µl	AT	4°C
AAC02-L	Acetyl-Lysine Antibody Mouse Monoclonal (7B5A1)	2 x 100 µl	AT	4°C
AAC02-beads	Acetyl-Lysine Affinity Beads - AAC02	4 x 500 ul	AT	4°C
AAC02-S	Acetyl-Lysine Antibody Mouse Monoclonal (7B5A1)	1 x 25 ul	AT	4°C
AAC03-L	Acetyl-Lysine Antibody Mouse Monoclonal (19C4B2-1)	2 x 100 µl	AT	4°C
AAC03-beads-L	Acetyl-Lysine Affinity Beads - AAC03	4 x 500 ul	AT	4°C
AAC03-HRP	Acetyl-Lysine-HRP Antibody Mouse Monoclonal (19C4B2-1)	1 x 100 ul	AT	4°C
AAC03-HRP-S	Acetyl-Lysine-HRP Antibody Mouse Monoclonal (19C4B2-1)	1 x 25 ul	AT	4°C
AAC03-S	Acetyl-Lysine Antibody Mouse Monoclonal (19C4B2-1)	1 x 25 ul	AT	4°C
AAC04-beads-A	Acetyl-Lysine Affinity Beads - AAC04	4 x 500 µl	AT	4°C
AACX1	Chemically Acetylated Bovine Serum Albumin (BSA)	2 x 500 µg	AT	4°C
AACX2	Tricostatin A (TSA) Treated And Untreated Cell Lysates	100 µg treated + 100	AT	4°C
APY03-L	Phosphotyrosine Antibody Mouse Monoclonal 27B10	2 x 100 µl	AT	4°C
APY03-Beads-L	Phosphotyrosine Affinity Beads	4 x 330 µl	AT	4°C
APY03-HRP-L	Phosphotyrosine-HRP Mouse Monoclonal Antibody 27B10	1 x 100 µl	AT	4°C
APY03-HRP-S	Phosphotyrosine-HRP Mouse Monoclonal Antibody 27B10	1 x 25 µl	AT	4°C
APY03-HRP-XL	Phosphotyrosine-HRP Mouse Monoclonal Antibody 27B10	Large Quantities	AT	4°C
APY03-S	Phosphotyrosine Antibody Mouse Monoclonal 27B10	1 x 25 µl	AT	4°C
ASM01-L	SUMO1 Antibody Mouse Monoclonal (Clone 5D8B16)	1 x 100 µl	AT	4°C
ASM01-S	SUMO1 Antibody Mouse Monoclonal (Clone 5D8B16)	1 x 25 µl	AT	4°C
ASM11-Beads-A	SUMOylation 1 Affinity Beads	3 x 400 µl	AT	4°C
ASM23-L	SUMO-2/3 Antibody Mouse Monoclonal (Clone 12F3)	2 x 100 µl	AT	4°C
ASM23-HRP-L	SUMO-2/3-HRP Antibody Mouse Monoclonal (Clone 12F3)	1 x 100 µl	AT	4°C

Cytoskeleton Inc. Campaign #4 2026
Signal Seekers

CAT. NO.	PRODUCT	AMOUNT	Shipping Temp.	Storage on Arrival
ASM23-HRP-S	SUMO-2/3-HRP Antibody Mouse Monoclonal (Clone 12F3)	1 x 25 µl	AT	4°C
ASM23-S	SUMO-2/3 Antibody Mouse Monoclonal (Clone 12F3)	1 x 25 µl	AT	4°C
ASM24-L	SUMO-2/3 Antibody Mouse Monoclonal (Clone 11G2)	2 x 200 µl	AT	4°C
ASM24-Beads-A	SUMOylation 2/3 Affinity Beads	2 x 400 µl	AT	4°C
ASM24-S	SUMO-2/3 Antibody Mouse Monoclonal (Clone 11G2)	1 x 150 µl	AT	4°C
AUB01-L	Ubiquitin Antibody Mouse Monoclonal	2 x 100 µl	AT	4°C
AUB01-HRP-L	Ubiquitin-HRP Antibody Mouse Monoclonal	1 x 100 µl	AT	4°C
AUB01-HRP-S	Ubiquitin-HRP Antibody Mouse Monoclonal	1 x 25 µl	AT	4°C
AUB01-S	Ubiquitin Antibody Mouse Monoclonal	1 x 25 µl	AT	4°C
AUB01-XL	Ubiquitin Antibody Mouse Monoclonal	4 x 500 µl	AT	4°C
BK160-L	Signal-Seeker® Phosphotyrosine Detection Kit	KIT 30 assays	AT	4°C
BK160-S	Signal-Seeker® Phosphotyrosine Detection Kit	KIT 10 assays	AT	4°C
BK161-L	Signal-Seeker® Ubiquitin Detection Kit	KIT 30 assays	AT	4°C
BK161-S	Signal-Seeker® Ubiquitin Detection Kit	KIT 10 assays	AT	4°C
BK162-L	Signal-Seeker® SUMOylation 2/3 Detection Kit	KIT 30 assays	AT	4°C
BK162-S	Signal-Seeker® SUMOylation 2/3 Detection Kit	KIT 10 assays	AT	4°C
BK163-L	Signal-Seeker® Acetyl-Lysine Detection Kit	KIT 30 assays	AT	4°C
BK163-S	Signal-Seeker® Acetyl-Lysine Detection Kit	KIT 10 assays	AT	4°C
BK164	Signal-Seeker® Acetyl-Lysine Detection Kit - AAC03-Beads (30 assays)	KIT 30 assays	AT	4°C
BK165-L	Signal-Seeker® SUMOylation 1 Detection Kit	KIT 30 assays	AT	4°C
BK165-S	Signal-Seeker® SUMOylation 1 Detection Kit	KIT 10 assays	AT	4°C
BLR01-50	BlastR Rapid Lysate Prep Kit	50 purifications	AT	AT
BLR02-50	BlastR Rapid Filtration Kit	50 purifications	AT	AT
BLR03-50	BlastR Lysis/Dilution Buffer Kit	50 lysates	AT	AT

Cytoskeleton Inc. Campaign #4 2026
Signal Seekers

CAT. NO.	PRODUCT	AMOUNT	Shipping Temp.	Storage on Arrival
CIG01-beads-10	IgG Control Beads	10 assays	AT	4°C
CIG02-beads-10	Acetyl-Lysine Control Beads	10 assays	AT	4°C
CIG03-beads-12	SUMO1 IP Control Beads	10 assays	AT	4°C
CUB02B-beads-A	Ubiquitination Control Beads	10 assays	AT	4°C
MB201-A	MsrB2 Protein 6xHis	2 x 50 µg	AT	4°C
MB201-XL	MsrB2 Protein 6xHis	1 x 1 mg	AT	4°C
MIC01-A	MICAL-1 Protein 6xHis	2 x 50 µg	AT	4°C
MIC01-XL	MICAL-1 Protein 6xHis	1 x 1 mg	AT	4°C
MXA95-A	MICAL-Oxidized Actin Protein (>95% pure) - Rabbit Skeletal Muscle	2 x 250 µg	AT	4°C
MXA95-XL	MICAL-Oxidized Actin Protein (>95% pure) - Rabbit Skeletal Muscle	1 x 1 mg	AT	4°C
MPAX1-A	MICAL-Oxidized (Pyrene labeled) Actin Protein (>95% pure) - Rabbit Skeletal Muscle	2 x 250 µg	AT	4°C
MPAX1-XL	MICAL-Oxidized (Pyrene labeled) Actin Protein (>95% pure) - Rabbit Skeletal Muscle	1 x 1 mg	AT	4°C
UBA01B-beads	Ubiquitination Affinity Beads	40 rxs	AT	4°C